

# Microzone Ltd

## *micro***CLEAN** for small DNA clean-up

**The easy way to clean-up or concentrate your PCR fragments and probes.**

### **Protocol:**

Just add equal volume of *micro***CLEAN** to DNA sample  
Mix by pipetting  
Leave at room temperature for 5 mins

### **For tubes:**


- Spin at high speed (10 to 13 K in microfuge) for 7 mins
- Remove supernatant
- Spin briefly again and remove all liquid (important)
- Resuspend pellet in the appropriate volume of TE
- Leave 5 mins to rehydrate DNA

### **For plates:**

- Spin at 2000 to 4000 g for 40 mins
- Place plate upside down onto tissue paper in the centrifuge holder
- Pulse centrifuge to <40 g for 30 secs
- Resuspend pellet in the appropriate volume of TE
- Leave 5 mins to allow DNA to rehydrate

**Clean DNA is now ready for further processing**

**Note: *micro*CLEAN-ed DNA is very clean and you should not see the pellet!**



Store at +4°C

For Research Only