


Microzone Ltd

100 bp DNA Ladder


Catalogue Number: 1DLS-100 and 1DLS-500
Volume: 1 x 1.5 ml for 100 loadings
5 x 1.5 ml for 500 loadings

100bp DNA
Ladder



The 100 bp DNA Ladder is supplied ready to use in loading buffer containing MZ Blue Dye. Bands of 80, 100, 200, 300, 400, 500, 600, 700, 800, 900, and 1000 bp

Load 10 to 15 μ l onto the gel (340 to 500 ng)



Store at +4 or -20°C

For Research Only