



# Microzone Ltd

## LYSE 'N EASY™

	<u>Volume</u>
microLYSISPLUS	1 ml
MegaMix-Royal	1 ml
Just water	1 ml

Enough for 50 x DNA release reactions  
and 100 x amplifications

Follow 2 stage protocol on a separate sheet

**(Note - microLYSIS-PLUS is a complex solution  
which releases (does not purify) the DNA from the  
cells. Estimation of DNA yield after lysis can  
therefore not be done using a spectrophotometer)**



Store at +4°C

For Research Only