


Microzone Ltd

Ethidium Bromide

10mg/ml
Molecular Biology Grade


Use at 0.5 $\mu\text{g/ml}$ for agarose gel electrophoresis
i.e. 5 μl stock solution per 100 ml agarose gel mixture or buffer.

Store at Room Temperature



Store at room temperature


For Research Only



Microzone Ltd

Mineral Oil

Molecular Biology Grade



Store at room temperature


For Research Only



MgCl₂


Molecular Biology Grade

MgCl₂ @ 1 M concentration



Store at room temperature


For Research Only



Microzone Ltd


Just Water

Molecular Biology Grade



Store at room temperature

For Research Only




Microzone Ltd

TE

**Molecular Biology Grade
10 mM Tris/1mM EDTA**

Store at Room Temperature



Store at room temperature

For Research Only