

# microzone


## 100 bp DNA Ladder

Ready to load agarose size standard

Containing 11 double stranded DNA fragments (34 ng/μl) in a room temperature stable loading buffer

### Description

Fragment sizes from top to bottom: 1000, 900, 800, 700, 600, **500**, 400, 300, 200, 100 & 80 bp



*A genuine Microzone product, made in the UK*

### Protocol

Load 10 to 15 μl per well in the agarose gel, alongside samples

### Related Products

**MegaMixes** = Clear or blue PCR mixes for instant and accurate PCRs

**1 Kb DNA ladder** = An agarose size standard supplied 'ready to load' in room temp. stable loading buffer

**microCLEAN** = The easy way to clean-up PCR fragments

Supplied as a blue solution in multiple of 1 ml aliquots  
Store at +4°C or at room temperature

For Research Only