

5x BaseQ[®] Sequence Dilution Buffer

Specially formulated to be used with ABI BigDye[®] Mix (versions 1.1 and 3.1)

Protocol

- Dilute the ABI BigDye Mix with BaseQ to achieve the desired concentration*
- Mix by pipetting
- Sequence according to protocol
- Remove excess dye

*1 in 2, 1 in 3, 1 in 4 or higher. The choice of dilution is dependent on size, quantity and quality of the DNA template. Please contact Microzone for answers to any questions
BigDye is a registered trademark of Applied Biosystems

Related Products

microCLEAN = The easy way to clean-up your PCR fragments

Just Water = Molecular grade water in handy sizes

Supplied as a clear solution in multiple of 1, 28 & 233 ml aliquots
Store at -20°C

For Research Only