

microFORCE[®]

Enhancing buffer for the amplification of GC rich DNA templates

Specially formulated to be used with MegaMix-Gold and MegaMix-Royal. Supplied at 4x concentration

Protocol

Use at 1x concentration (i.e. 5 µl microFORCE in 20 µl final volume PCR)

Related Products

MegaMix-Gold and MegaMix-Royal = Double concentrated HotStart PCR mixes

100bp DNA ladder = An agarose size standard supplied 'ready to load' in room temperature stable loading buffer

Just Water = Molecular grade water in handy sizes

Supplied as a clear solution in multiple of 1 ml aliquots
Store at +4°C

For Research Only